



Massage chair A2801













Armrest shortcut key
USB charging



Speaker
Home sound-transmitting cloth

Product Configuration:

Apperance: PU+PVC

Rail: SL, Zero space, dual electric cylinder

Massage robot: 2D

Footrest: open toe, spring telescopic

Control system: convenient remote controller+voice control+armrest shortcut key

Heating: waist

Charging: A/C port

Speaker: bluetooth

Air pressure: Airbag pressure detection, adjusts the inflation time according to body shape automatically.

Safety: backrest anti-pinch detection


Brackets: remote controller bracket

Packing: Disassemble the calf and put it on the seat cushion and pack it as one.

60 units/40HQ



Leisure massage chair
A2801



Airbag pressure detection, adjusts the inflation time according to body shape automatically.

Continuously monitors the air pressure within the airbags. When the pressure reaches the preset value, the air valve automatically ceases inflation. The shoulder, arm, leg, and foot airbags feature multi-level pressure adjustment, better accommodating the needs of different individuals.




Zero Gravity


Position the thighs and calves horizontally above the heart, allowing the body's weight to distribute evenly across the massage chair and achieve a zero-gravity state. With the body's centre of gravity resting on the hips, pressure on the spine and joints is reduced, enabling complete physical relaxation.



Accurate shoulder position detection

After the massage chair is automatically lying in place, it will automatically detect the human shoulder position, and the acupuncture point positioning will be more accurate, and the massage robot will perform precise massage according to the body shape to get the best massage physiotherapy effect.






The enlarged heating area for the waist, combined with kneading massage for the lower back, accelerates blood circulation in the lumbar region to effectively alleviate back pain. It achieves constant-temperature heat therapy at 40°C–50°C within three minutes. During autumn and winter, the heat application function promotes blood circulation in the feet, offering beneficial health effects for women, the elderly, and those prone to feeling cold, providing users with warm care.

Bluetooth music

Immersive stereo surround music: The left and right sides of the massage chair hood are equipped with a stereo Bluetooth speaker, surround music, enjoy music while massaging, create a comfortable massage atmosphere, so that you can relax your whole body.





Upright: 140.5*75.5*119.5cm

Lying down: 166.5*76*94.5cm

Packing: 1 CTN

Package dimension: 116.5*77.5*126cm

Loading quantity: 60pcs/HQ

Rated voltage: 110V-240V

Rated power: 80W

Calf lift angle: 95°

Maxmium lying down angle: 155°

Distance from wall: 5cm

N.W.: 约85Kg

G.W.: : 约100Kg

THANKS!